

AS and A-level Chemistry practicals: Equipment set up

Practical 4: Carry out simple test-tube reactions to identify

- cations Group 2, NH⁴+
- anions Group 7 (halide ions), OH–, CO₃²–, SO₄²–

Test 3: testing for ammonium ions

- The tube should be placed in a warm water bath and a piece of damp red litmus held near the mouth of the tube using tongs.
- The red litmus should change colour and become blue in the presence of ammonium ions.





Test 4: hydroxide ions

Damp red litmus paper changes to blue.



Ammonium ions in a petri dish.

The damp red litmus is pressed on the inside of the lid.





Pipette collection of carbon dioxide for test 5

Several pipettes of carbon dioxide will be needed to show a change in the colour of the limewater.





Collection of carbon dioxide for test 5 using a delivery tube


- The delivery tube can be glass or plastic and fastened in the neck of the test tube with a rubber bung.
- Limewater in a second tube is used to test for carbon dioxide.





Halide tests


Results of the test for halides with the addition of 2M ammonia solution.





Halide tests


Results of the test for halides after the addition of concentrated ammonia solution.





Halide tests

Results of the halides with silver nitrate solution.





Results of halides with concentrated sulfuric acid.

This must be done in a fume cupboard.



Potassium bromide result





Potassium chloride result





Potassium iodide result



